

Frankel Appliances



Frankel I Appliance: Corrects overcrowding in Class I cases and reduces the overbite and overjet in Class II, Division 1 cases. External muscle pressure is eliminated by the vestibular shields. Promotes transverse arch development dentally and skeletally.

Frankel II Appliance: Promotes transverse and vertical development of maxillary and mandibular arches. Corrects Class II, Division 2 cases and opens bite. Used after the maxillary incisors have been slightly proclined by an upper removable appliance.

Frankel III Appliance: Corrects Class III malocclusions.

Frankel IV Appliance: Used for open bites and bimaxillary protrusions.

Address: 5/F, Building 3, Goldstar Industrial Park, Jianan Road, Fuhai Street, Baoan District, Shenzhen China

Tel: +86-13828833659

E-mail: info@wmdentallab.com

website: www.wmdentallab.com